Please Note

There is a newer version of this document available at:

RelianceOrthodontics.com/ifu

RelianceOrthodontics.com

RelianceOrthodontics.com

RelianceOrthodontics.com

RelianceOrthodontics.com

RelianceOrthodontics.com

RelianceOrthodontics.com

RelianceOrthodontics.com

RelianceOrthodontics.com

RelianceOrthodontics.com

RelianceOrthodontics.com

RelianceOrthodontics.com

RelianceOrthodontics.com

RelianceOrthodontics.com

RelianceOrthodontics.com

RelianceOrthodontics.com

RelianceOrthodontics.com

RelianceOrthodontics.com

RelianceOrthodontics.com

RelianceOrthodontics.com

RelianceOrthodontics.com

RelianceOrthodontics.com

RelianceOrthodontics.com

RelianceOrthodontics.com

RelianceOrthodontics.com

RelianceOrthodontics.com

RelianceOrthodontics.com

RelianceOrthodontics.com

RelianceOrthodontics.com

RelianceOrthodontics.com

RelianceOrthodontics.com

RelianceOrthodontics.com

RelianceOrthodontics.com

RelianceOrthodontics.com

RelianceOrthodontics.com

RelianceOrthodontics.com

RelianceOrthodontics.com

RelianceOrthodontics.com

RelianceOrthodontics.com

RelianceOrthodontics.com

RelianceOrthodontics.com

RelianceOrthodontics.com

RelianceOrthodontics.com

RelianceOrthodontics.com

RelianceOrthodontics.com

RelianceOrthodontics.com

RelianceOrthodontics.com

RelianceOrthodontics.com

EN - ENGLISH

INDICATIONS FOR USE: Maximum Cure® Sealants are intended for use in direct and indirect bonding of brackets and appliances.

INSTRUCTION FOR USE:

CUSTOM BASE INDIRECT BONDING

- Prophy the tooth surface. Rinse and dry. Acid etch each tooth for 30 seconds. Rinse and dry.
- Apply 1 coat of , Assure® Plus to teeth and each custom pad, air dry.
- There are 2 choices for use as an indirect transfer adhesive:
 - Mix Maximum Cure part A & B and apply a thin coat to the enamel. Then, apply a slightly thicker coat to each custom pad and immediately seat tray.
 - OR
 - Apply a moderate coat of Maximum Cure part A to the custom pad and part B to the enamel. Proceed to seat the tray.

In either case, hold the tray for 90 seconds and wait 5 minutes before removing tray from mouth.

Special Safety Notice: Brush applicators are single use devices. It is required that a new brush applicator is used for each application to prevent contamination.

⚠️⚠️ **PRECAUTIONS:** Unpolymerized resins may cause skin sensation in susceptible persons. In case of contact with the skin, wash thoroughly with soap and water. If accidentally instilled into eye, flush with copious amounts of water and seek medical attention immediately. Contact with other tissues: rinse spill area immediately with copious amounts of water for several minutes.

⚠️⚠️ **Flammable** ⚠️ **Hazardous to aquatic life**

⚠️⚠️ **Corrosive**

🔥🔥 **TEMPERATURE LIMITATION:** Store refrigerated (2-8°C).

Avoid excessive heat.

TO REORDER: #MC7AF, #MC7B, #MC20AF, #MC20B, #MC77AF, #MCF5B

WARRANTY: Reliance Orthodontics Products, Inc. recognizes its responsibility to replace products if proven to be defective. Reliance Orthodontics Products, Inc. does not accept liability for any damages or loss, either direct or consequential, stemming from the use of or inability to use the products as described. Before using, it is the responsibility of the user to determine the suitability of the product for its intended use. The user assumes all risk and liability in connection therewith. **RX ONLY:** U.S. Federal law restricts this device to sale by or on the order of a dental professional.

NOTE: Report any serious incident to Reliance and the competent authority of the state in which the user and/or patient is established.

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