

ReliaBox Quick Start Guide



# Unlock the benefits of ReliaBox

The key to efficient pre-pasted bracket placement



**ReliaBox innovative storage solution is designed to revolutionize your bracket application process.**

Pre-paste brackets for patient appointments up to a week in advance, reducing chairside time and increasing convenience for you and your team.

### ✔ Smart & Systematic Design.

ReliaBox is divided into 4 quadrants with 4 anterior tweezers and 3 posterior tweezers per quadrant.

### ✔ Accessible & Versatile.

When hung correctly, tweezers protrude 1-2" for convenient chairside access and hangers are designed to be used with any bracket holders, not just the Reliance brand.

### ✔ Protects Against Ambient Air & Light.

The custom designed locking seal allows for clinicians to pre-paste brackets up to 7 days in advance - keeping out both ambient air and light.

### ✔ Efficient & Effective Chairside.

Shave at least 15 minutes off every chairside appointment utilizing the ReliaBox pre-paste system!  
\*IMPORTANT: Do not place used tweezers back in box. Place with other used instruments for sterilization.

See it in action

