




1540 West Thorndale Ave., Itasca, IL 60143 U.S.A.
1-800-323-4348 · www.RelianceOrthodontics.com

**SELF MIXING DUAL-SYRINGE DELIVERY SYSTEM
INSTRUCTIONS FOR USE**

1. Remove cap from the SELF MIXING SYRINGE by rotating in a counter-clockwise direction and pulling off.
2. Place mixing tip by lining up notch in tip hub (A) with the syringe notch (B) and rotating in a clockwise direction until it is in a locked position.
3. For first time use of SELF MIXING SYRINGE: Dispense the first 2-3mm of material onto a mixing pad to ensure a homogeneous mixture of material.



4. Syringe material into band/crown or appliance using mixing tip.
5. **DO NOT PULL BACK ON THE PLUNGER AFTER DISPENSING.** This can cause cross contamination of the base and catalyst, curing the paste and plugging the syringe.
Do not re-use the tips. The SELF MIXING SYRINGE mixing tip should be removed and discarded after each patient use. If used intraorally, the syringe should be cold sterilized prior to storage.
6. Replace the cap on the syringe after use.
CAUTION: DO NOT REUSE IF SYRINGE BECOMES CONTAMINATED. AVOID CONTACT OF THE REUSABLE PARTS WITH THE PATIENTS MOUTH. USE NEW, UNCONTAMINATED GLOVES. APPLY DISPOSABLE BARRIER SLEEVES/WRAPS BEFORE EACH PATIENT


WARRANTY: Reliance Orthodontic Products, Inc. recognizes its responsibility to replace products if proven to be defective. **Reliance Orthodontic Products, Inc.** does not accept liability for any damages or loss, either direct or consequential, stemming from the use of or inability to use the products as described. Before using, it is the responsibility of the user to determine the suitability of the product for its intended use. The user assumes all risk and liability in connection therewith.



1540 West Thorndale Ave., Itasca, IL 60143 U.S.A.
1-800-323-4348 · www.RelianceOrthodontics.com

**SELF MIXING DUAL-SYRINGE DELIVERY SYSTEM
INSTRUCTIONS FOR USE**

1. Remove cap from the SELF MIXING SYRINGE by rotating in a counter-clockwise direction and pulling off.
2. Place mixing tip by lining up notch in tip hub (A) with the syringe notch (B) and rotating in a clockwise direction until it is in a locked position.
3. For first time use of SELF MIXING SYRINGE: Dispense the first 2-3mm of material onto a mixing pad to ensure a homogeneous mixture of material.



4. Syringe material into band/crown or appliance using mixing tip.
5. **DO NOT PULL BACK ON THE PLUNGER AFTER DISPENSING.** This can cause cross contamination of the base and catalyst and plugging of the syringe.
Do not re-use the tips. The SELF MIXING SYRINGE mixing tip should be removed and discarded after each patient use. If used intraorally, the syringe should be cold sterilized prior to storage.
6. Replace the cap on the syringe after use.
CAUTION: DO NOT REUSE IF SYRINGE BECOMES CONTAMINATED. AVOID CONTACT OF THE REUSABLE PARTS WITH THE PATIENTS MOUTH. USE NEW, UNCONTAMINATED GLOVES. APPLY DISPOSABLE BARRIER SLEEVES/WRAPS BEFORE EACH PATIENT

WARRANTY: Reliance Orthodontic Products, Inc. recognizes its responsibility to replace products if proven to be defective. **Reliance Orthodontic Products, Inc.** does not accept liability for any damages or loss, either direct or consequential, stemming from the use of or inability to use the products as described. Before using, it is the responsibility of the user to determine the suitability of the product for its intended use. The user assumes all risk and liability in connection therewith.